



## EZ-iPSC Generation Kit - EPISOMAL

### Product Information

**Catalog Number** ASK-3013 (Episomal vector based iPSC reprogramming kit)

**Description** Induced pluripotent stem cells (iPSCs) are genetically reprogrammed from adult cells, which are similar to natural pluripotent stem cells, such as embryonic stem cells (ESCs). iPSCs exhibit a pluripotent stem cell-like state, such as the expression of certain stem cell genes and proteins, chromatin methylation patterns, doubling time, teratoma formation, and potency and differentiability. While these artificially generated cells are not known to exist in the human body, they show qualities remarkably similar to those of embryonic stem cells. Therefore, iPSCs are an invaluable resource for drug discovery, cell therapy, and basic research.

The Human EZ-iPSC Generation Episomal Kit is an optimized mixture of multiple vectors that can reprogram somatic cells to iPSCs without integration. The episomal vectors have the oriP/EBNA-1 (Epstein-Barr nuclear antigen-1) backbone that delivers the reprogramming factors (Oct-3/4, Sox2, Klf4, and l-Myc) as well as puromycin resistance gene. This system has been successfully demonstrated in the reprogramming of fibroblasts, as well as other adult cells. High expression of transgenes due to oriP/EBNA-1 mediated nuclear import and retention of vector DNA allows iPSC derivation from a single transfection. The reprogramming efficiency is further enhanced by puromycin selection. In addition, silencing of the viral promoter which drives EBNA-1 expression and the loss of episomes at a rate of ~5% per cell cycle due to partitioning and defects in vector synthesis allows the removal of episomal vectors from the iPSCs without any additional manipulation.

The episomal vectors are introduced into the cell by electroporation and the reprogramming protocol is quite simple yet highly effective with efficiencies in the range of 0.05-1.0%, depending upon the cell type being reprogrammed. Reprogramming enhancers, such as PD0325901, CHIR99021, A-83-01, hLIF, and HA-100, are not required for reprogramming.

**Contents** Each kit contains sufficient material for about 10 reprogramming experiments:  
EZ-iPSC Generation episomal vectors (ASK-3013-1): 2 vials, 25  $\mu$ L each (1  $\mu$ g/ $\mu$ L)  
RFP Control vector (ASK-3013-2): 1 vial, 10  $\mu$ L (0.5  $\mu$ g/ $\mu$ L)

**Shipping** Blue ice

**Storage** Store at -20 °C upon receipt.

**Quality Control** Each lot of EZ-iPSC Generation Episomal Kit is tested to ensure human fibroblasts are efficiently reprogrammed into iPSCs (figure 1), for absence of mycoplasma and viral pathogens, HIV 1, HIV II, Hepatitis B, and Hepatitis C.

**Biosafety Level** BSL-2 laboratory

Applied StemCell, Inc.

521 Cottonwood Dr. #111, Milpitas, CA 95035

Phone: 866-497-4180 (US Toll Free); 408-773-8007 Fax: 408-773-8238

[info@appliedstemcell.com](mailto:info@appliedstemcell.com) [www.appliedstemcell.com](http://www.appliedstemcell.com)


Copyright 2018, Applied StemCell, Inc. All rights reserved. This information is subject to change without notice.

**Safety Precaution**

**PLEASE READ BEFORE HANDLING ANY FROZEN VIALS.** Please wear the appropriate Personal Protection Equipment (lab coat, thermal gloves, safety goggles and a face shield) when handling the cells. Handle the frozen vials with due caution. Please be aware that the following scenario can occur: Liquid nitrogen can leak into the vials when the vials are submerged in liquid nitrogen. Upon thawing, the liquid nitrogen returns to the gas phase, resulting in a dangerous build-up of pressure within the vial. This can result in the vial exploding and expelling not only the vial contents but also the vial cap and plastic fragments of the vial.

**Restricted Use**

This product is for research use only and not intended for human or animal diagnostic or therapeutic uses.



**Figure 1.** Immunostaining of iPSC colonies reprogrammed from human fibroblasts using the Human EZ-iPSC Generation Episomal kit (ASK-3013): iPSC colonies expressing stem cell markers specific to Oct4 and TRA-1-60.

**Media and Material**

EZ-iPSC Generation Episomal Kit (ASK-3013) contents: sufficient for 10 reprogramming experiments.

Catalog #	Component	Product Details	Storage
ASK-3013-1	EZ-iPSC Generation Episomal vectors	2 vials: each 25 $\mu$ L (1 $\mu$ g/ $\mu$ L)	-20°C
ASK-3013-2	RFP control vector	1 vial: 10 $\mu$ L (0.5 $\mu$ g/ $\mu$ L)	-20°C

**Materials and Instruments required but not provided**

Materials Needed	Vendor	Quantity
Neon transfection devices and kits	Invitrogen	100 $\mu$ L sets
iPSGen Medium	Stem RD	
PSGro Medium	Stem RD	
mTeSR1	Stem Cell Technologies	
Matrigel	BD	

## Protocol

For optimal reprogramming with the Human EZ-iPSC Generation Episomal Kit, culture the fibroblasts in Fibroblast Medium until the day of transfection. After transfection, allow the cells to recover in Fibroblast Medium for 24 hours, then add puromycin to remove the untransfected cells. Reseed the cells on feeders about 5 day of post transfection, and switch the culture medium to hESC culture medium the following day.

### 1. Feeder-free reprogramming protocol

#### 2.1 Preparation of Fibroblast Medium with pen/strep

DMEM containing 10% FBS, 2mM GlutaMax, and 50 U and 50  $\mu$ g/ ml penicillin and streptomycin.

#### 2.2 Preparation of Fibroblast Medium without pen/strep

DMEM containing 10% FBS, 2mM GlutaMax.

#### 2.3 Preparation of PSGro/mTeSR1 Medium

1:1 mix PSGro medium (StemRD) and mTeSR1 (Stemcell Technologies), and add penicillin and streptomycin to final concentration of 50 U and 50  $\mu$ g/ ml, respectively.

#### 2.4 Coating plates with Matrigel

Matrigel should be aliquoted and stored at -80°C for long-term use.

1.4.1 Thaw Matrigel on ice until liquid. Dilute Matrigel 1:50 with pre-chilled KO DMEM/F12.

1.4.2 Immediately use the diluted Matrigel solution to coat tissue culture-treated plates. For a 6-well plate, use 0.8 ml of diluted Matrigel solution per well, and swirl the plate to spread the Matrigel solution evenly across the surface.

1.4.3 Let the coated plate stand for 1 - 2 hrs at 37°C or overnight at 4°C. If plate has been stored at 4°C, allow the plate to incubate at 37°C for at least 30 minutes before removing the Matrigel solution.

#### 2.5 Procedure



Day -1: Seed Cells

- 1.5.1 Seed human fibroblasts at  $1 \times 10^6$  cells in a gelatin coated T75 flask. Cells should reach approximately 75–90% confluent on the day of transfection (Day 0).

## Day 0: Nucleofection of Episomal Plasmids to Human Fibroblasts

- 1.5.2 Prepare the Neon Transfection Devices and kits (using 100- $\mu$ l tips in this protocol).
- 1.5.3 Prepare the gelatin-coated 6 well plate and warm up Fibroblast Medium (w/o P/S) in the plate.
- 1.5.4 Aspirate the medium from the fibroblasts in the T75 flask, and wash the cells in DPBS without calcium and magnesium. Add 2 ml of 0.05% Trypsin/EDTA to the flask, and incubate the flask at 37°C for 3 minutes.
- 1.5.5 Add 5 ml of Fibroblast Medium to the flask. Tap the flask to ensure the cells are dislodged from the flask, and carefully transfer the cells into a new 15-ml conical tube. Spin down the cells at 1,000 rpm for 5 minutes at room temperature.
- 1.5.6 Resuspend the cells in 2 ml of DPBS. Count the cells, and take  $3 \times 10^5$  cells into each 1.5 ml tube for two tubes. Note: You need  $3 \times 10^5$  cells for one transfection, and the cell number can be modified to 1.5 – 6  $\times 10^5$  cells per transfection.
- 1.5.7 Spin down the cells at 2000 rpm for 5 min. In the meantime, add 3ml of solution E2 to the microporation tube. Mix 5  $\mu$ g of reprogramming vectors in one 1.5 ml tube with 100  $\mu$ l Solution R. In another tube, mix 1.5  $\mu$ g for RFP control with 100  $\mu$ l Solution R. Note: If you use 10- $\mu$ l tip, mix 1.5  $\mu$ g of reprogramming vectors in 10  $\mu$ l Solution R and 0.5  $\mu$ g for RFP control in 10  $\mu$ l Solution R, respectively.
- 1.5.8 Carefully aspirate most of the supernatant. Resuspend the cell pellet in Solution R with plasmids.
- 1.5.9 Turn on the Neon unit and enter the electroporation parameters in the Input window to 1,650 Volts of pulse voltage, 10 ms of pulse width, and 3 pulses. Note: To increase the viability, you may use the electroporation parameters of 1,400 Volts of pulse voltage, 30 ms of pulse width, and 1 pulse.
- 1.5.10 Press the push-button on the Neon® Pipette to the first stop and immerse the Neon® Tip into the cell-DNA mixture. Slowly release the push-button on the pipette to aspirate the cell-DNA mixture into the Neon® Tip. Note: Avoid air bubbles when pipetting to avoid arcing during electroporation. If you notice air bubbles in the tip, discard the sample and carefully aspirate fresh sample into the tip again without any air bubbles.
- 1.5.11 Insert the Neon® Pipette with the sample vertically into the Neon® Tube placed in the Neon® Pipette Station until you hear a click. Ensure that you have entered the appropriate electroporation parameters and press Start on the Neon® touchscreen to deliver the electric pulse.
- 1.5.12 Immediately after electroporation, the cell suspension solution is poured into warm Fibroblast Medium w/o P/S) in one well of a 6-well plate pre-coated with gelatin. Culture the electroporated cells in 37 °C, 5% CO2 incubator.

## Day 1, 3, 5: Change the Fibroblast Medium

- 1.5.13 Change medium to Fibroblast Medium with P/S, supplemented with 0.5  $\mu$ g/ml of puromycin.

*Note: If RFP plasmid is used to monitor the transfection efficiency, the RFP expression can be detected during these days. By adding puromycin, non-puro-resistant cells dislodge.*

*Note: Do NOT treat transfected cells with puromycin for more than three days.*

## Day 6: Replating the Transfected Cells

- 1.5.14 Aspirate the medium from the well with transfected fibroblasts, and wash the cells in DPBS without calcium and magnesium. Add 0.5 ml of 0.05% Trypsin/EDTA to the well, and incubate the plate at 37°C for 2 minutes.
- 1.5.15 Add 1 ml of Fibroblast Medium to the well. Tap the plate to ensure cells are dislodged from the plate, and carefully transfer cells to a new 15-ml conical tube. Spin down cells at 1,000 rpm for 5 minutes at room temperature.
- 1.5.16 Resuspend the cells in 2 of ml Fibroblast Medium. Count the cells. Seed  $1 \times 10^4$ ,  $2 \times 10^4$ , and  $4 \times 10^4$  transfected fibroblasts, respectively, in the wells of a 6-well plate pre-coated Matrigel. Note: The cell numbers can be changed from 0.5 -  $3 \times 10^5$  cells when using a 100-mm dish. The seeding day can vary according to the cells.

**Day 7: Switch to iPSGen Medium**

- 1.5.17 Aspirate the medium from the reprogramming fibroblasts in the well, and add 2 ml of iPSGen Medium (StemRD). Note: You may use N2B27 medium supplemented with 100 ng/ml of bFGF to replace iPSGen Medium.
- 1.5.18 Change the medium everyday up to day 14.

**Day 14: Switch to PSGro/mTeSR1 Medium**

- 1.5.19 Aspirate the medium from the reprogrammed fibroblasts, and add 2 ml of PSGro/mTeSR1 medium. Note: You may use mTeSR1 medium to replace PSGro/mTeSR1 medium.
- 1.5.20 Change PSGro/mTeSR1 medium every day.

**Day 24 – 30: Picking iPS-like Colonies**

By Day 24 of post-transfection, the cell colonies in the Matrigel-coated plate consist of iPSCs, which exhibit a hESC like morphology characterized by a flatter, cobblestone-like appearance with individual cells clearly demarcated from each other in the colonies.

- 1.5.21 Examine the culture plate containing the reprogrammed cells under 10X magnification of an inverted microscope, and mark the colony to be picked at the bottom of the culture dish. Note: We recommend picking at least 10 distinct colonies by the end of each reprogramming experiment and expanding them in separate 24-well Matrigel-coated plate.
- 1.5.22 Transfer the culture dish to a biosafety cabinet equipped with a stereomicroscope.
- 1.5.23 Cut the colony to be picked into 5–6 pieces in a grid-like pattern using a 25-gauge 1½ inch needle.
- 1.5.24 Using a 200-µl pipette, transfer the cut pieces to a Matrigel-coated well of a 24-well plate containing mTeSR1 medium supplemented with 10 µM ROCK inhibitor (Y-27632, StemRD).
- 1.5.25 Culture the picked colonies in a 37°C, 5% CO<sub>2</sub> incubator.
- 1.5.26 Allow the colonies to attach to the culture plate for 48 hours before replacing the spent medium with fresh mTeSR1 medium. Since then, change the medium every day.
- 1.5.27 Culture the reprogrammed colonies like normal human iPSC colonies; expand and maintain them using standard culture procedures.

*Note: Newly derived iPSC lines may contain a fair amount of differentiation through passage 4. It is not necessary to remove differentiated material prior to passaging. By propagating the cells, the overall culture health should improve throughout the early passages. Otherwise, pick up the iPS-like colonies and culture in Matrigel coated 24-well plate.*

**2. Reprogramming Protocol Using MEF Feeders**

Materials Needed: Human iPS Cell Reprogramming Episomal Vectors 25 µl (1 µg/µl)  
 Episomal RFP control vector 10 µl (0.5 µg/µl)  
 Neon Transfection Devices and kits (Invitrogen; 100 µl sets)  
 Feeder cells (mitomycin-c treated-SNL or MEF cells)

**2.1 Preparation of Fibroblast Medium with pen/strep**

DMEM containing 10% FBS, 2mM GlutaMax, and 50 U and 50 µg/ ml penicillin and streptomycin.

**2.2 Preparation of Fibroblast Medium without pen/strep**

DMEM containing 10% FBS, 2mM GlutaMax.

**2.3 Preparation of human ES Medium**

Knockout DMEM/F12 containing 20% knockout serum replacement, 2mM glutamine, 0.1 mM nonessential amino acids, 0.1 mM 2-mercaptoethanol, 10 ng/ml bFGF, and 50 U and 50 µg/ ml penicillin and streptomycin.

**2.4 Preparation of hES/PSGro Medium**


1:1 mix PSGro Medium (StemRD) and human ESC Medium, and add penicillin and streptomycin to final concentration of 50 U and 50 µg/ ml, respectively.

## 2.5 Preparation of mitomycin-C treated MEF or SNL feeders

MEF cells can be treated with mitomycin C when they reach confluency, to halt the division of cells.

- 2.5.1 Add 6 ml of fresh MEF (or SNL) medium contain 50  $\mu$ l of mitomycin C solution (1 mg/ml) to one T75 flask of confluent MEF (or SNL) cells, and swirl it briefly. The final concentration of mitomycin C is 8  $\mu$ g/ml.
- 2.5.2 Incubate at 37°C for at least 3 hrs.
- 2.5.3 Aspirate the mitomycin C-containing medium off the cells and wash the cells twice with 10 ml PBS.
- 2.5.4 Aspirate PBS and add 2 ml of 0.25% trypsin-EDTA, swirl to cover the entire surface, and incubate for 2 min at room temperature.
- 2.5.5 Add 5 ml Fibroblast Medium and break up the cells to single-cell suspension by pipetting up and down. Count the number of cells.
- 2.5.6 Seed the cells on gelatin-coated dishes (3 x 10<sup>6</sup> cells per 100-mm dish, or 5 x 10<sup>5</sup> cells per well of a 6-well plate).
- 2.5.7 Cells should be ready to use on the next day.

## 2.6 Procedure



### Day -1: Seed Cells

- 2.6.1 Seed human fibroblasts at 1 x 10<sup>6</sup> cells in a gelatin coated T75 flask. Cells should reach approximately 75–90% confluence on the day of transfection (Day 0).

### Day 0: Nucleofection of Episomal Plasmids to Human Fibroblasts

- 2.6.2 Prepare the Neon Transfection Devices and kits.
- 2.6.3 Prepare the gelatin-coated 6 well plates and warm up Fibroblast Medium (w/o P/S) in the plates.
- 2.6.4 Aspirate the medium from the fibroblasts in the T75 flask, and wash the cells with DPBS without calcium and magnesium. Add 2 ml of 0.05% Trypsin/EDTA to the flask, and incubate the flask at 37°C for 3 minutes.
- 2.6.5 Add 5 ml of Fibroblast Medium to the flask. Tap the flask to ensure the cells are dislodged from the flask, and carefully transfer the cells into a new 15-ml conical tube. Spin down the cells at 1,000 rpm for 5 minutes at room temperature.
- 2.6.6 Resuspend the cells in 2 ml of DPBS. Count the cells, and take 3 x 10<sup>5</sup> cells into each 1.5 ml tube for two tubes. Note: You need 3 x 10<sup>5</sup> cells for one transfection, and the cell number can be modified to 1.5 – 6 x 10<sup>5</sup> cells per transfection.
- 2.6.7 Spin down the cells at 2000 rpm for 5 minutes. In the meantime, add 3ml of solution E2 to the microporation tube. Mix 3 $\mu$ g of reprogramming vectors in one 1.5 ml tube with 100  $\mu$ l Solution R. In another tube, mix 1  $\mu$ g for RFP control with 100  $\mu$ l Solution R. Note: If you use 10- $\mu$ l tip, mix 1.5  $\mu$ g of reprogramming vectors in 100  $\mu$ l Solution R and 0.5  $\mu$ g for RFP control in 10  $\mu$ l Solution R, respectively.
- 2.6.8 Carefully aspirate most of the supernatant, and resuspend the cell pellet in Solution R with plasmids.
- 2.6.9 Turn on the Neon unit and enter the electroporation parameters in the Input window to 1,650 Volts of pulse voltage, 10 ms of pulse width, and 3 pulses. Note: To increase the viability, you may use the electroporation parameters of 1,400 Volts of pulse voltage, 30 ms of pulse width, and 1 pulse.
- 2.6.10 Press the push-button on the Neon® Pipette to the first stop and immerse the Neon® Tip into the cell-DNA mixture. Slowly release the push-button on the pipette to aspirate the cell-DNA mixture into the Neon® Tip. Note: Avoid air bubbles when pipetting to avoid arcing during electroporation. If you



notice air bubbles in the tip, discard the sample and carefully aspirate fresh sample into the tip again without any air bubbles.

- 2.6.11 Insert the Neon® Pipette with the sample vertically into the Neon® Tube placed in the Neon® Pipette Station until you hear a click. Ensure that you have entered the appropriate electroporation parameters and press Start on the Neon® touchscreen to deliver the electric pulse.
- 2.6.12 Immediately after electroporation, the cell suspension solution is poured into warm Fibroblast Medium (w/o P/S) in a 6-well plate pre-coated with gelatin. Culture the electroporated cells in 37 °C, 5% CO<sub>2</sub> incubator.

## Day 1, 3, 5: Change the Fibroblast Medium

- 2.6.13 Change medium to Fibroblast Medium with P/S, supplemented with 0.5 µg/ml of puromycin. Note: If RFP plasmid is used to monitor the transfection efficiency, the RFP expression can be detected during these days. By adding puromycin, non-puro-resistant cells dislodge.
- 2.6.14 On day 5, seed the MEF/SNL cells in a gelatin-coated plate or dish (5 x 10<sup>5</sup> cells per well of a 6-well plate or 3 x 10<sup>6</sup> cells per 100-mm dish).
- 2.6.15 Day 6: Replating the Transfected Cells
- 2.6.16 Aspirate the medium from the transfected fibroblasts in the well, and wash the cells in DPBS without calcium and magnesium. Add 0.5 ml of 0.05% Trypsin/EDTA to the well, and incubate the plate at 37°C for 2 minutes.
- 2.6.17 Add 1 ml of Fibroblast Medium to the well. Tap the plate to ensure the cells have been dislodged from the plate, and carefully transfer the cells into an empty 15-ml conical tube. Spin down the cells at 1,000 rpm for 5 minutes at room temperature.
- 2.6.18 Resuspend the cells in 2 ml of Fibroblast Medium. Count the cells. Seed 1 x 10<sup>4</sup>, 2 x 10<sup>4</sup>, and 4 x 10<sup>4</sup> transfected fibroblasts, respectively, onto wells of a 6-well plate pre-seeded MEF/SNL feeder cells. Note: The cell numbers can be changed from 5 x 10<sup>4</sup> to 3 x 10<sup>5</sup> cells when using a 100-mm dish. The seeding day can vary according to the cell types.

## Day 7: Switch to hES Medium

- 2.6.19 Aspirate the medium from the fibroblasts in the well, and add 2 ml of hES Medium supplemented with 10 ng/ml of bFGF.
- 2.6.20 Change the medium everyday up to day 15.

## Day 16: Switch to hES/PSGro Medium

- 2.6.21 Aspirate the medium from the fibroblasts in the well, and add 2 ml of hES/PSGro Medium (StemRD). *Note: You may also use mTeSR1 (Stemcell Technologies) or Essential 8 Medium (Life Technologies) to replace PSGro medium.*
- 2.6.22 Change hES/PSGro Medium every day.

## Day 21 – 30: Picking iPSC-like Colonies

By Day 21 of post-transfection, the cell colonies in the MEF plate consist of iPSCs, which exhibit a hESC-like morphology characterized by a flatter, cobblestone-like appearance with individual cells clearly demarcated from each other in the colonies (see Figure 2).

- 2.6.23 Examine the culture plate containing the reprogrammed cells under 10X magnification of an inverted microscope, and mark the colony to be picked at the bottom of the culture dish. Note: We recommend picking at least 10 distinct colonies by the end of each reprogramming experiment and expanding them in separate 24-well Matrigel-coated plates.
- 2.6.24 Transfer the culture plate to a biosafety cabinet equipped with a stereomicroscope.
- 2.6.25 Cut the colony to be picked into 5–6 pieces in a grid-like pattern using a 25-gauge 1½ inch needle.
- 2.6.26 Using a 200-µL pipette, transfer the cut pieces to a well of a 24-well plate pre-seeded MEF/SNL feeder cells. Note: If you want to culture the iPSCs in feeder-free conditions, transfer the cut pieces to a Matrigel-coated well of a 24-well plate containing mTeSR1 medium supplemented with 10 µM ROCK inhibitor (Y-27632, StemRD).
- 2.6.27 Incubate the picked colonies in a 37°C, 5% CO<sub>2</sub> incubator.

- 2.6.28 Allow the colonies to attach to the culture plate for 48 hours before replacing the spent medium with fresh hES medium. After that, change the medium every day.
- 2.6.29 Treat the reprogrammed colonies like normal human iPSC colonies; expand and maintain them using standard culture procedures.

*Note: Newly derived iPSC lines may contain a fair amount of differentiation through passage 4. It is not necessary to remove differentiated material prior to passaging. By propagating the cells, the overall culture health should improve throughout the early passages. Otherwise, pick the iPS-like colonies and culture in 24-well plate with feeder cells.*