



Jurkat-C9m (Cas9 expressing Jurkat Cells)

Order information


Catalog Number	Quantity
ASE-9300	5.0 x 10 ⁵ cells/ vial

Product Information

Description

Guide RNA (gRNA) validation is usually performed in K562 or HEK293 cell lines. The gRNA activity validated in these cell lines cannot always be translated in blood cells. A Jurkat cell line was engineered to stably express the *S. pyogenes* Cas9 nuclease. This Cas9 expressing Jurkat cell line can be used for gRNA validation by transfecting the gRNA into the cell line followed by mismatch assay.

The data below indicate the validation result of the guide RNAs. The Jurkat-C9m was transfected with three oligo-synthesized guide RNAs. Genomic DNA was harvested 72 hours after transfection. All three guide RNAs are shown to be active in the mismatch assay.



Shipping

Dry ice

Storage and Stability

Store in liquid nitrogen freezer immediately upon receipt. This product is stable for at least 6 months from the date of receiving when stored as directed.

Quality Control


To ensure the expression of Cas9, the cell line should be periodically cultured in hygromycin 100ug/ml. The positive guide RNA (gPC) can be used as a control for gRNA validation assay.

Applied StemCell, Inc.

521 Cottonwood Dr. #111, Milpitas, CA 95035

Phone: 866-497-4180 (US Toll Free); 408-773-8007 Fax: 408-773-8238

info@appliedstemcell.com www.appliedstemcell.com



Safety Precaution

PLEASE READ BEFORE HANDLING ANY FROZEN VIALS. Please wear the appropriate Personal Protection Equipment (lab coat, thermal gloves, safety goggles and a face shield) when handling the cells. Handle the frozen vials with due caution. Please be aware that the following scenario can occur: Liquid nitrogen can leak into the vials when the vials are submerged in liquid nitrogen. Upon thawing, the liquid nitrogen returns to the gas phase, resulting in a dangerous build-up of pressure within the vial. This can result in the vial exploding and expelling not only the vial contents but also the vial cap and plastic fragments of the vial.

Restricted Use

This product is for research use only and not intended for human or animal diagnostic or therapeutic uses.

Protocol

1. Cell culture

Media Component	Final Conc.
RPMI 1640 (supplemented with FCS)	
Fetal calf serum (FCS)	10%
Streptomycin	100 µg/mL
Penicillin	100 U/mL
GlutaMAX™	2 mM

- 2.1 Seeding density: 1 x 10⁵ cells/ mL
 2.2 Passage interval: Cells should be passaged before reaching density of 4-6 million cells/mL

2. Guide RNA preparation

- 2.1 Oligo-synthesis: Order synthetic guide RNA
 2.2 *In vitro* transcription: Use HiScribe T7 High Yield RNA Synthesis Kit (NEB, Cat. #: E2040S) to transcribe gRNAs.
- 2.2.1 Clean the lab bench and pipetors by wiping with RNaseZap® RNase Decontamination Solution.
 - 2.2.2 Keep the T7 RNA Polymerase Mix on ice.
 - 2.2.3 Thaw the T7 10X Reaction Buffer, four ribonucleotide solutions, and ddH₂O at room temperature.
 - Briefly vortex the T7 10X Reaction Buffer and ribonucleotide solutions. Centrifuge before open tubes.
 - 2.2.4 Prepare template:
 - HiScribe T7 High Yield RNA Synthesis Kit suggests 1 µg template DNA for short RNA transcription (0.3 kb).
 - Ambion MEGAShortscript Kit suggests 260 ng DNA templates in 20 µL reaction to yield best amount of RNA (~90 µg) from 100 bp-length DNA.

- 2.2.5 Assemble the reaction mix in an RNase-free microfuge tube at room temperature as shown below.
 2.2.6 dNTPs, 10X Buffer and Enzyme Mix can be premixed as 9 μ L mixture added into 20 μ L reaction):
- For short RNA transcripts (0.3 kb)

Segments	Volume (μ L)
T7 10X Reaction Buffer	1.5
dNTPs mix (75 mM)	1.5 X 4
T7 Enzyme Mix	1.5
Template DNAs (100 ng - 1 μ g)	11

- 2.2.7 Mix contents thoroughly by gently flicking the tube.
 2.2.8 Incubate the reaction at 37°C for at least 4 hrs.

3. Transfection: Lipofectamine

- 3.1 Plate 2-3 x 10⁵ cells/ml density of cells with fresh medium in one well of 12-well plate the day before transfection.
 3.2 Add 2.5 μ L of Dharmacon gRNA plus 2.5 μ L of Dharmacon rRNA to 50 μ L OPTI-MEM and add 3 μ L Lipofectamin 3000 to 50 μ L OPTI-MEM. Mix together for 5 min.
 3.3 Add 2 μ g IVT gRNA to 50 μ L OPTI-MEM and add 3 μ L Lipofectamin 3000 to 50 μ L OPTI-MEM. Mix together for 5 min.
 3.4 Add 100 μ L transfection solution to cells
 3.5 Collect the transfected cells 72 hours after transfection for mismatch essay. Refresh the medium if necessary during the culture.

4. Mismatch Assay

- 4.1 Extract genomic DNA.
 4.2 PCR 100ng of genomic DNA using 2x Fusion (Fisher, Cat #: F548L)
 4.3 Run 2 μ L of PCR sample on 2% agarose gel with Low Mass ladder (Invitrogen, Cat. #:10068-013) to estimate concentration.
 4.4 Adjust concentration to 400ng in 17 μ L volume with diH₂O.
 4.5 Hybridize PCR sample to form heteroduplexes using program below in a thermocycler.

Temperature	Time	Temperature ramp
95°C	10 min	
95°C to 85°C		(-2.0°C/s)
85°C	1 min	
85°C to 75°C		(-0.3°C/s)
75°C	1 min	
75°C to 65°C		(-0.3°C/s)
65°C	1 min	
65°C to 55°C		(-0.3°C/s)
55°C	1 min	
55°C to 45°C		(-0.3°C/s)
45°C	1 min	
45°C to 35°C		(-0.3°C/s)
35°C	1 min	
35°C to 25°C		(-0.3°C/s)
25°C	1 min	
4°C	Hold ∞	

- 4.6 Add 2 μ L of 10X NEBuffer 2 and 1 μ L of T7 Endonuclease I (NEB, cat. #M0302S)

4.7 Incubate at 37°C for 15min. Add 2 μ L of 0.5M EDTA to stop the reaction.

4.8 Run 16 μ L of sample on acrylamide gel.

5. Cryopreservation

Resuspend the cells in the culture medium with supplement of 10% DMSO and place vials into an isopropanol freezing container and place the container at -80°C overnight. Transfer to a liquid nitrogen tank on the following day.