



## Stromal Cells Culture Protocols

### Product Information

Catalog Number **ASM-6004 (Stromal Cellulations Medium)**

### Materials Required

Materials	Equipment
<ul style="list-style-type: none"><li>• Cells ASE-5005, ASE-5006, ASE-5007, ASE-5194^5196</li><li>• 1X PBS (Ca<sup>2+</sup>/Mg<sup>2+</sup>-free)</li><li>• Stromal Cellulations Media (ASM-6004) Collagen coating solution (150µg/mL)</li><li>• Dimethyl Sulfoxide (DMSO)</li><li>• Fetal Bovine Serum (FBS)</li><li>• TrypLE Express</li><li>• 15.0 &amp; 50.0 mL tubes</li><li>• Tissue culture vessels or 6 well plates</li><li>• 2-1000µl Pipets</li><li>• 2-100µl Pipet Tips</li><li>• Pipet Aid</li><li>• Serological Pipet Tips</li><li>• Mr. Frosty</li><li>• Ice</li></ul>	<ul style="list-style-type: none"><li>• Incubator, 37°C/5% CO<sub>2</sub></li><li>• Water Bath 37°C</li><li>• Centrifuge</li><li>• Pipette</li></ul>

### Protocol

#### 1. Immediately Upon Delivery

- 1.1 Remove vial from shipping container to check that it is still frozen.
- 1.2 Transfer frozel vial to liquid nitrogen until you are ready to thaw and begin cell culture.

#### 2. Thawing Cells

- 2.1 Prepare complete Stromal Cellulations Media (can be stored at 4°C for a maximum of 30 days).
- 2.2 Aliquot and warm only media required.
- 2.3 Thaw cells rapidly and with agitation (See Hot to Thaw Cells).
- 2.4 Centrifuge at 400 g for 5 minutes. Remove supernatant and resuspend cell pellet in 1-5 mL fresh 37°C warm Stromal Cellulations Media.
- 2.5 Count cells (see How to Count Fresh Cells).

#### 3. Cultivating Cells

- 3.1 Prepare an aliquot of cells for plating at a density:
  - 1.0-5.0 x 10<sup>4</sup> cells/cm<sup>2</sup> for mononuclear cells
  - 1.0 x 10<sup>3</sup> cells/cm<sup>2</sup> for stromal cells.
- 3.2 Gently mix warm media with cells for plating.
- 3.3 Carefully add cell suspension prepared to flasks/plate at correct seeding density.
- 3.4 Place flasks/plates in a CO<sub>2</sub> incubator.
  - For mononuclear cells, let cells sit undisturbed for no less than 3-4 days, change media on the 5<sup>th</sup> day.
  - For Stromal Cells, proceed to next step.

#### Applied StemCell, Inc.

521 Cottonwood Dr. #111, Milpitas, CA 95035

Phone: 866-497-4180 (US Toll Free); 408-773-8007 Fax: 408-773-8238

[info@appliedstemcell.com](mailto:info@appliedstemcell.com) [www.appliedstemcell.com](http://www.appliedstemcell.com)

*Note: You will lose a lot of non-adherent mononuclear cells at this point. See note below for continuing with non-adherent cells.*

- Feed cells every 2-3 days by performing a complete media change as follows:
- Aspirate Media
  - Add warm PBS
  - Aspirate PBS
  - Add fresh media

*Note: Mononuclear cells will have the vast majority of the cells floating. They will not adhere. Collect the supernatant if you are interested in investigating other cell types besides stromal.*

- Visualize cells every 2-3 days.
- Once cells become 70% confluent prepare cells for passage, collection, or freezing.

*Note: Mononuclear cells will take 3-4 days to attach and may take up to 3 weeks to reach 70% confluent.*

*Note: Do not let mononuclear cells grow to more than 70% confluency or for > 14 days.*

#### 4. Passing and Sub-cultivation of Cells

4.1 Aspirate and discard medium; wash cells with PBS (Ca<sup>2+</sup>/Mg<sup>2+</sup> free); add TrypLE enzyme solution or similar and incubate for 5-10 min at 37°C. (for details see How to Passage Cells).

*Note: Visualize cells and collect as soon as they are lifted off the plate. DO NOT over trypsinize.*

4.2 Transfer all enzyme solution mixture into a conical tube containing 10% volume of cold H-GRO or FBS, keep it chilled on ice.

4.3 Wash flasks/plate with equal volume PBS and collect any remainder cells

4.4 Centrifuge the collected cell mixture in the conical tube at 400 g for 5 minutes, and aspirate supernatant

4.5 If plating cells to another cell culture vessel:

- Resuspend cell pellet in pre-warmed media
- Perform cell count
- Plate cells at  $1.0 \times 10^3$  cells/cm<sup>2</sup>
- Refer to product information sheet for specifics on characterization of cells.

4.6 If freezing cells for storage:

- Resuspend cell pellet in cold FBS
- Remove aliquot and perform cell count
- Place cells on ice
- Determine how many vials to freeze (recommended  $1.0 \times 10^6$  to no more than  $5.0 \times 10^6$  in 1mL per vial)
- Slowly and dropwise add equal volume cold FBS containing 20% (v/v) dimethyl sulfoxide (DMSO). This achieves a final volume of 10% (v/v) dimethyl sulfoxide (DMSO).
- Immediately remove vials to Mr. Frosty and place at -80°C or control rate freezer.