


iPSC Generation, Disease Modeling, Mutation Correction and Differentiation

Comprehensive suite of premium quality stem cell related services and products for affordable and efficient modeling of “disease-in-a-dish”

Why work with ASC?

- ASC is one of the earliest licensees and most experienced service provider of both the iPSC (from iPS Academia, Japan) and the CRISPR/Cas9 technology (from the Broad Institute), worldwide.
- ASC’s CRISPR and iPSC disease modeling platforms have been featured in top peer-reviewed journals such as Nature, and PNAS, as well as in biotechnology magazines.
- ASC’s genome editing services include engineering cell lines from healthy and diseased donors.
- iPSC genome editing portfolio includes customizable deliverables, including a full panel of pluripotency characterization and teratoma analysis.



Footprint-Free iPSC Generation Service

ASC has successfully reprogrammed hundreds of iPSCs from both healthy and disease patient samples, using its proprietary footprint-free, non-integrating protocols. Some key features include:

- Reprogramming from patient specific fibroblasts or blood cells
- No genomic footprint: suitable for drug discovery and cell therapy applications
- High reprogramming efficiency
- Safe to handle: no viral particles

Case Study: Generation of human iPSCs from dermal fibroblasts of a patient with a rare disorder. The dermal fibroblasts obtained by skin biopsy from a patient were reprogrammed into iPSC using ASC’s proprietary protocols. High quality colonies were isolated and characterized by immunocytochemical staining for pluripotency markers.




Figure. iPSC reprogrammed from patient-derived fibroblasts and cultured using ASC’s proprietary protocols. The reprogrammed iPSC colonies stained positive for pluripotency markers OCT4, SOX2, SSEA4, TRA-1-60, and TRA-1-81, and alkaline phosphatase (AP).

You need to provide

- 1 x 10⁶ cells of patient fibroblast or PBMC
- Pathogen-free report

You will receive

- At least 3 clones with 2 vials of each clone (> 2 x 10⁵ cells/vial)
- Full report with high resolution images
- Clones are characterized by pluripotency markers
- Teratoma analysis, EB formation, qPCR, RNA-seq (optional)

Timeline: 2 months

Cell recovery and quarantine



Transfection with episomal DNA vectors



Colony formation



Colony expansion



Characterization via antibody staining and RT-PCR (Teratoma and EB analysis optional)




Cell Banking/further differentiation



iPSC Genome Editing/Disease Modeling Service (Knock-in, Knockout, Point Mutation)

ASC's acclaimed iPSC gene editing service provides precise engineering of iPSCs to your specifications, using highly optimized iPSC and CRISPR/Cas9 transfection conditions to achieve >90% success rate. With our unique expertise in genome engineering and stem cells, our custom iPSC engineering offers the following advantages:

- Isogenic control and disease cell lines for reliable comparison of results without genetic variability
- Unlimited resource of *in vitro* models of human genetics and diseases
- Patient-related research and drug discovery



iPSC Genetic Modification

- CRISPR/Cas9
- TARGATT™

ASC provides services under licensing agreement with the Broad Institute for CRISPR/Cas9 Gene-Editing technology and iPS Academia, Japan for iPSC technology along with ASC's proprietary TARGATT™

You need to provide

- 1 x 10⁶ cells of your iPSC or using ASC's iPSC lines
- Gene information

You will receive

- Clones with desired mutation
- Full report with publishable data

- Customized Deliverables:**
- Choice of homozygous and/or heterozygous clones
 - Choice to use silent mutations
 - Choice to use feeder-free or feeder-dependent protocols
 - Foot-print free editing

- Gene knockout, insertion, and replacement
- Reporter gene insertion
- Gene correction using your patient-derived iPSCs, or ASC's characterized iPSCs

Case Studies

1. Knockout Mutation

Goal: To generate a frameshift knockout (KO) mutation in ASC's human iPSC cell line using CRISPR/Cas9. ASC's iPSC line was used for this study. After gRNA validation and transfection, the iPSCs were screened and potential KO iPSC colonies were subjected to genotyping PCR and sequencing.




Figure 1. Genotyping analysis and sequence confirmation of KO iPSC clones. PCR amplification of either side of the gRNA cut site indicated the presence of a 551 bp PCR product with a 22 bp deletion clone #6; wildtype (WT) PCR product = 573 bp. (B) Sequence chromatogram showing deletion. (C) Sequence alignment between and wildtype.



iPSC Disease Cell Line Model Generation Service

Based on our proprietary site-specific TARGATT™ technology, ASC developed a human TARGATT™ iPSC line with pre-engineered attB sites for introducing gene(s) of interest into human iPS cells. The unique advantage of this system is that any gene of interest can be inserted efficiently into a defined, transcriptionally-active locus with high gene expression in a single copy fashion.

Start from iPSC with TARGATT™ site:

TARGATT™ iPSCs allow easy insertion of any gene of interest into a safe harbor locus in the genome.

Site-specific rapid knockin:

- High knock-in efficiency
- Single **transgene** copy
- Minimal off-target

Timeline: 3-5 months

Cell recovery and quarantine

TARGATT DNA vectors construction and validation

Transfection of TARGATT constructs

Cell confirmation and expansion





Figure 2. Staining of TARGATT™ iPSC master cell line with standard iPSC expression markers. (a) Oct4, (b) Sox2, (c) SSEA4, (d) TRA-1- 60, (e) TRA-1-81.

You need to provide

- Gene of interest
- Project details

You will receive

- Clones with desired gene knock-in
- Full report with publishable results

TARGATT™ iPSC Master Cell Line & Quick Knock-in Kit

ASC's TARGATT™ iPSC Quick Knock-in Kit is a complementary kit that works with our TARGATT™ iPSC Master Cell Line to insert large fragment DNA into preselected "attP" docking sites that have been engineered into the iPSC Master Cell Line. This kit is an excellent tool for generating cell line models for gene overexpression, and to generate reporter lines in hiPS and differentiated cells. The TARGATT™ iPSC Genotyping Kit provides a convenient method to genotype the TARGATT™ iPSC Master Cell Line or Knock-in cells generated from a Master Cell Line.

- AST-1100 iPSC Master Cell Line
- AST-1101 Quick knockin kit
- AST-1102 Genotyping kit now available

Knock-in Kit Contains

- Non-replicable PhiC31 integrase plasmid
- Positive control plasmid with GFP
- TARGATT™ cloning plasmids (with different promoters)

Genotyping Kit Contains

- DNA extraction buffers
- Two sets of genotyping PCR primers flanking integration sites and inserted fragment
- Nuclease free water
- Positive control DNA



iPSC Characterization Services

Teratoma Formation Analysis Service

Teratoma Analysis is one of the most stringent and accurate quality control assessments for stem cells, in addition to other methods to characterize stemness such as immunohistochemistry and RT-PCR detection of pluripotency markers, in vitro embryoid body formation and karyotyping. It provides a functional assessment of pluripotency of the stem cells, including iPSCs, by analyzing the ability of the cells to form all three embryonic germ cell layers when xenografted into mice. ASC's teratoma formation analysis service has >97% success rate and has been acknowledged in > 30 peer-reviewed publications.



You will need to provide

- Two T75 flasks containing cells for testing (note: total count of 3-6 x 10⁶ cells for mouse iPSCs).

Timeline: 1-3 months

Cell injection
(2 sites: kidney and testis)

Teratoma harvesting

Tissue sectioning

H&E staining

Histology analysis of teratoma sections




Figure 1. Histological analyses using H&E staining of kidney and testis teratomas from mice injected with an iPSC line. Differentiated tissues representing the three germ layers are shown and indicated by arrow heads. Abbreviations: EN: endoderm; ME: mesoderm; EC: ectoderm. Magnification: 200x.

You will receive

- Complete report with histological analysis and high quality/ publication grade images of EN, ME and EC formation.
- Tissue blocks and H&E stained tissue section slides.

Success Rate

Key points	ASC's methods			Traditional methods
Cell Type	mESC/miPSC	hESC	hiPSC	hESC
Teratoma Formation Rate	100%	100%	93.7%	25-40%
Differentiation	distinctive	distinctive	distinctive	poor
Turnaround Time	3-5 weeks	5-8 weeks	10-14 weeks	12-18 weeks
Cells needed	0.5-1 million/site	0.5-2 million/site	1-2 million/site	3-5 million/site

Antibody Pluripotency Marker Staining:

- Human: OCT4, SOX2, SSEA4, TRA-1-60, TRA-1-81
- Mouse: OCT4, SOX2, SSEA1

Pluripotency & lineage-specific marker detection qPCR, RNA-seq

Karyotyping: Chromosome counting; G-banding)

Embryoid Body Formation and Characterization

Germline Transmission Evaluation

mESC Derivation Service

Custom iPSC Culture Service

iPSC Banking Service

Related Services* (*additional fees may apply)

Cell Immortalization Service

- Patient fibroblasts, primary cells, and more
- Transduction, colony selection (up to 10 passages)
- Characterization (transgene expression by RT-PCR)


Vector Cloning, Virus Packaging Services

- RNAi and inducible vectors/ BAC recombineering
- Retro- and lentivirus packaging service in 10 days (regular, high, ultra-high titers)



iPSC Differentiation Service

ASC offers various iPSC differentiation services for ESC / iPSC differentiation into more specialized cells including multipotent stem cells and fully differentiated somatic cells. Using proprietary protocols and reagents, we can differentiate your ESCs or patient-derived or healthy iPSCs into various lineage cells: neural stem cells (NSC), neurons, cardiomyocytes, and hepatocytes



- You need to provide**
- Cells to be differentiated: 2 vials of frozen stocks on dry ice or 1 T25 flask of live cells per sample
 - Pathogen-free report

- You will receive**
- Fresh or frozen cells
 - Comprehensive characterization of differentiated cells
 - Full report with high resolution images

Human Neural Stem Cell Differentiation Service

ASC will use proprietary protocols and optimized reagents to differentiate of your iPSCs into ready-to-use multipotent neural stem cells. Simply thaw the cells we ship and you can use them directly for your neuroscience research or further differentiate into motor neurons and/or oligodendrocytes.

- High purity (>90%) NSCs in 1 week using ASC's proprietary neural induction protocol
- Derived NSCs are proliferating cells that can be frozen and cultured over a prolonged period of time while retaining phenotype
- Derived NSCs retain phenotype under feeder-free culture conditions
- NSCs can be further differentiated into multiple neural cell types.




Figure. iPSC differentiation service. (A) Time line to show NSC differentiation workflow. (B) iPSCs differentiation to NSCs in feeder-free condition with high efficiency using ASC's proprietary neural induction protocol and Neural Differentiation media (ASM-4013) in < 1 week. The differentiated NSCs retain neural stem cell phenotype even after passaging and freeze-thaw under feeder-free conditions as evidenced by staining for NSC markers PAX6 and SOX2 and nuclear DNA (DAPI). Note: White arrows indicate distinctive neuronal rosettes.

NSC-related products: Neural Stem Cells and Medium

- ASE-9303 Human iPSC-derived Neural Stem Cells
- ASM-4013 NeuroSure™ Neural Differentiation Medium
- ASM-4014 NeuroSure™ Neural Stem Cell Culture Medium



Human Cardiomyocyte Differentiation Service

ASC provides custom service to differentiate iPSCs into functional cardiomyocytes from your iPSCs directly, or after re-programming your patient-derived fibroblasts and PBMCs (optional custom service). Our proprietary protocols yield high purity (>90%), beating cardiomyocytes that express typical markers such as TNNT/cTnT and α -Actinin (see Figure 1), and are validated for functional viability by patch-clamp and Fluo-4TM Direct Calcium Assay. These cardiomyocytes can be used for electrophysiological and biochemical assays (Figure 2). For example, the effect of compounds on the heart rate can be assayed by impedance-based measurements.




Figure 1. After replating in Cardiomyocyte Maintenance Medium on the Matrigel-coated plates for 2 days, more than 85% of the iPSC-derived cardiomyocytes express cTnT and α -Actinin (Figure A and B).




Figure 2. Beating rate of iPSC-derived cardiomyocytes measured with xCELLigence RTCA Cardio System (ACEA Biosciences and Roche Applied Science). Cells were treated with various dosages of Isoproterenol. (Left) Before Treatment (50 μ M; 10 μ M; 5 μ M); (Right) 30 min after treatment (50 μ M; 10 μ M; 5 μ M).

iPSC Cells from ALS8 Patients and Non-Carrier Siblings

ASC's iPSC catalog now includes fully characterized iPSCs reprogrammed from fibroblasts of 5 individual ALS8 patients (with P56S mutation) and their non-affected siblings belonging to two families affected with the autosomal dominant form of familial ALS (ALS8). This series of iPSC lines is a very useful and new tool to study ALS with reduced genetic variability and can be used to develop early diagnostic tools, and to identify possible targets for drugs and therapies.





Figure. Characterization of controls and ALS8-patients iPSCs. (A) Heredograms of the two ALS8 families. ALS8 patients: dark symbols; non-affected individuals: white; *: Skin biopsies from affected and control individuals. (J-M) In vivo characterization by teratoma induction identified tissue from the 3 germ layers.

Ref: Mitne-Neto, et al. (2011). Downregulation of VAPB expression in motor neurons derived from induced pluripotent stem cells of ALS8 patients. *Human molecular genetics*, 20(18), 3642-3652.

iPSC Cell Culture Products

MEF Cells

- Strict quality control tested for long term culture of either mouse or human ESC/iPSC lines (>168 passages)
- Our MEFs are used by more than 200 clients in academic and industrial labs worldwide
- DR4, CF-1, Neo, SNL 76/7 cells available (Untreated, irradiated, mitomycin C-treated)

Cell Culture Products

- Serum and feeder free medium
- Conditioned medium for human ESC/iPSC culture
- Freezing medium

Easy-to-use Kits


- EZ iPSC generation kit (Episomal/Retroviral)
- iPSC characterization kit
- Immortalization kit

Primary Cells

- Fibroblasts (Healthy/Disease)
- Cardiomyocytes, Cardiac Progenitor Cells

ESC-Sure™ FBS

- High quality, affordable ESC-grade serum for ES and iPS cell culture
- Supports undifferentiated growth of mESCs (QC report provided)



ESC-Sure™ FBS is extensively tested for supporting undifferentiated growth of mouse ESCs (mouse embryonic stem cells). The image shows healthy mESC that were cultures in MEF-feeder cell dependent system along with ESC's ES-grade FBS.

We provide high-quality MEF cells!
Average viability >95%



My EZ Gel™ CDX Models

- Protein-based solubilized base membrane
- Synthesized from commonly used amino acids
- Compatible with various adhesion proteins and cell growth factors
- Powerful reagent to package and deliver cancer cells, primary cancer cells, or stem cells

ESC/ iPSC Cell Line Catalog

Human iPSC Cells from Healthy and Patient Samples

- ASE-9203 Human iPSCs (Episomal, Healthy, Fibroblasts)
- ASE-9028 Human iPSCs (Episomal, Parkinson's Disease, Fibroblasts)
- ASE-9029 Human iPSCs (Episomal, Parkinson's Disease, Adipocytes)
- ASE-9030 Human iPSCs (Episomal, Parkinson's Disease, PBMCs)
- ASE-9031 Human iPSCs (Episomal, Diabetes II, Fibroblasts)
- ASE-9033 Human iPSCs (Episomal, Diabetes II, PBMCs)

iPSCs & ESCs: Mouse, Rat, Pig, Guinea Pig

- ASE-9005 Germline mESC; C57Bl/6
- ASE-9006 Germline mESC; 129-EZ
- ASE-9007 Germline mESC C57Bl/6 EZ
- ASE-9008 Germline mESC BALB/c EZ
- ASE-9009 mESC C57Bl/6 Albino
- ASE-9106 Pig iPSCs
- ASE-9107 Mouse iPSCs
- ASE-9108 Guinea Pig iPSCs

