



INTRODUCTION


- CRISPR gene editing in cell lines, the workhorse of preclinical and biomedical research, enables the generation of unlimited *in vitro* models with precise gene modifications and advanced gene expression design that are physiologically relevant.
- CRISPR genome editing in certain types of cells, such as suspension cells and several blood lineage cell lines is very inefficient and problematic.
- Contributing factors include low transfection efficiency of CRISPR reagents, cytotoxicity, and low or undesired gRNA and Cas9 activities in the cell lines.
- In this poster, we describe techniques optimized for CRISPR genome editing in some difficult-to-transfect blood-lineage cell lines, the factors influencing efficiency such as cytotoxicity, and the types of modifications achieved (double gene knockout, large fragment knock-in, point mutation).
- These results demonstrate that different cell lines may require different approaches or modified protocols to deliver CRISPR components for efficient and successful modifications of the targeted gene.

CRISPR/Cas9 Cell Line Engineering



Large Fragment Deletion in Jurkat Cell Lines


Using Cas9-plasmid:



(A) Schematic representation of ~5kb targeted deletion for enhancer region of the GOI (gene of interest) generated by co-delivery of two gRNAs (g1 and g2). The primer design for detection of region of deletion is indicated by orange arrows. (B) Two single clones were examined by PCR and gel electrophoresis, followed by Sanger sequencing. Clone #2 showed a shorter band at the size of 300bp and was confirmed to be positive for the deletion of the ~5 kb fragment in the gene of interest. (C) Sequence alignment of target region between positive clone (#2) and wild type (WT) indicates deletion of required sequence. DNA ladder: GeneRuler 1kb plus DNA ladder.

Efficient HDR-Mediated Homozygous Point Mutations in K562 Cells


Using a DNA-dependent protein kinase inhibitor to improve efficiency of HDR mediated sequence insertion:



(A) Schematic representation of the donor sequence (REF-seq) that will serve as the HDR template to introduce the CAT > AAC point mutation in the gene of interest in K562 cell lines. The gRNA, PAM and silent mutation targeting designs are represented accordingly. A silent mutation was introduced to prevent re-cleavage by Cas9. (B) CRISPR/Cas9-mediated HDR efficiency is improved upon co-treatment with a DNA-dependent protein kinase inhibitor for 24h. (C) Sequence alignment of representative clones indicate required point mutation in the targeted region. Treatment with the drug generated indels in both alleles of the targeted region in the cell line. Representative heterozygous and homozygous clones are also aligned with reference sequence to indicate required mutation in the targeted locus. WT: wild type sequence; REF-seq: reference sequence/donor sequence with silent mutation.


Gene Knockout in Leukemia Cell Lines, KG-1 Cells

Using Cas9-RNP:




Double Gene Knockout in Jurkat Cell Lines

Using Cas9-plasmid:



Knock-in & Knockout in T2 Lymphoblastic Cell Line

Using DNA-based CRISPR reagents:



(A) Schematic representation of CRISPR/Cas9 targeting strategy for transgene (gene of interest; GOI) knock-in into a safe harbor locus in Lymphoblastic cell line, T2. (B) PCR gel electrophoresis of single cell clones indicate clones A1, A2, A3 are positive for knock-in of GOI at targeted ROSA26 locus as compared to wild type (WT). PCR with flanking primers to confirm insertion at ROSA26 locus (a) and (b) show fragment of required sizes: Left arm (LA): ~2.4 kb, Right arm (RA): ~2.0kb. PCR to confirm transgene knock-in, MA also shows fragment of required size: ~2.2kb. Lane M: GeneRuler 1kb plus DNA ladder.

CONCLUSION

- CRISPR/Cas9 is a powerful tool for engineering many different cell lines from several species.
- CRISPR/Cas9 enables the generation of knockout, knock-in and point mutation lines with precise genetic modifications.
- CRISPR efficiency can be limited by the type or cell line being engineered, especially difficult-to-transfect cell lines such as blood-lineage cell.
- A variety of factors such as cell growth parameters, cell line sensitivity and cytotoxicity to transfected CRISPR reagents, can impact efficiency and reliability in genome editing.
- Optimizing protocols and transfection conditions for these blood lineage cells can contribute to higher cell survival, efficiency and success rate.
- Our results demonstrate that difference cell lines may require different approaches or modified protocols to deliver CRISPR/Cas9 components for efficiency and success modifications of targeted genes.