



Combined CRISPR - TARGATT™ Gene Editing Generates High Efficiency Knock-out, Point Mutation and Large Fragment Knock-in in Human iPS Cells

Luping Huang, Andreia Sommer, Padmaja Tummala, Lingjie Kong and Ruby Yanru Chen-Tsai

Applied StemCell Inc., Milpitas, CA 95035

Introduction


Gene editing in iPS cells offers great promise for personalized therapy in regenerative medicine

- CRISPR works most efficiently for gene knockout through the non-homologous end joining (NHEJ) pathway, and for point mutation/correction through homology directed repair (HDR).
- CRISPR efficiency for large fragment DNA knock-in through nuclease-mediated HDR is very low.
- Applied StemCell's TARGATT™ integrase-based system can successfully knock-in large fragment DNA in iPS cells (up to 22 kb)
 - Site specific gene insertion with high efficiency (up to 40%)
 - Transgene integration happens at a pre-selected, transcriptionally active safe-harbor locus
- **Combined CRISPR - TARGATT™ technology broadens scope and efficiency of gene editing in human iPS cells**


Materials and Methods

- Human iPS cells (ASC Cat# ASE-9202)
- CRISPRCLEAR™ hROSA26 Safe Harbor Knock-in Kit (ASC Cat# ASK-7041)
- Electroporation using Life Technology Neon® Transfection System

CRISPR & TARGATT™ Genome Editing Technologies




TARGATT™ Integrase System Site-specific Gene Knock-in in iPS Cells



Site-Specific Knock-in Technology
TARGATT™


CRISPRCLEAR™

Knockout Mutation in Six Isoforms of a Gene in Human iPS Cells



Six Isoforms of the gene of interest shared the same ATG start codon in targeted exon. (A) gRNA target site in gene of interest. (B) Deep sequencing was used to identify the best gRNA candidate, g2 with lowest off-target and highest on-target efficiency. (C) Co-transfection of gRNA and Cas9 in hiPSCs resulted in an **8 bp deletion** in one allele and a **1 bp insertion** in second allele at the DSB site, resulting in premature stop codons in all six isoforms of the gene of interest.

Point Mutation Correction in Human Disease iPS Cell Line




Homozygous point mutation correction in a disease human iPS Cell Line. (A) gRNA selection strategy in gene of interest. (B) Sequence chromatogram confirming the **GTA > ATA point mutation correction** in the disease cell line. The selected gRNA was transfected into patient cell line along with an ssODN containing the corrected wild type codon. Shown is a homozygous clone with the corrected mutation.

Combining CRISPR & TARGATT™ to Enable Large DNA Insertion

TARGATT™ Master Cell Line: Using CRISPR to Insert a "docking site"



TARGATT™ Master Cell Line: Large DNA Insertion



Conclusions

- Our data show that CRISPR can enable to generation of knockout and point mutation human iPS cell line models.
- Applied StemCell's TARGATT™ technology enables efficient knock-in of large DNA fragments in iPS cells at the preselected Rosa26 locus.
- Combining CRISPR & TARGATT™ technologies provides a valuable new platform for efficient genome editing in human iPS cells.