



# Custom FFPE Block Order Form

For Order Assistance Please Contact Technical Support at:

Tel.1-866-497-4180 (US Toll Free) Tel.1-408-773-8007

Fax.1-650-800-7179 email form to: [orders@oncoref.com](mailto:orders@oncoref.com)

for quote: save pdf with new file name and click button

1) **Customer Info:** Name \_\_\_\_\_ email \_\_\_\_\_  
Institution \_\_\_\_\_ Telephone \_\_\_\_\_

2) **Choose Cell Line for preparing FFPE Cell Block**

a. Do you want to use your own cell line? Yes or No If No, skip to Section 1b.

If yes, will you need us to grow the cell line for you Yes or No

b. Do you want to use [ONCOREF™ Mutation Isogenic Cell Line](#)? Yes or No

If yes, list ONCOREF™ Catalog Number \_\_\_\_\_ (e.g. ASE-8023)

c. Do you want to use Human Cancer Cell Lines?

If yes, choose from list below \_\_\_\_\_ (e.g. SK-BR-3)

- Breast Cancer (SK-BR-3, MCF7, MDA-MB-231, MDA-MB-468, BT-474, T-47D)
- Lung Cancer (H2228, H23, A549, H460, HOP-62, H322, SK-MES-1)
- Cervical Cancer (HeLa, Ca Ski)
- Prostate Cancer (PC3, DU-145)
- Colon Cancer (SW-620, HT-29, HCC-2998)
- Melanoma (SK-MEL-2, SK-MEL-5, M14, M21, CHL-1)
- Leukemia (K-562, Jurkat, Raji)
- Ovarian Cancer (OVCAR-8, IGR-OV1)
- Renal Cancer (SN-12C, TK-10)
- Liver Cancer (HepG2, Huh-7)

2) **Custom FFPE block preparation**

a. How many FFPE Blocks need to be prepared? \_\_\_\_\_

b. Do you need multiple (i.e duplicate) blocks from the same cell line? Yes  or No .

If yes, how many? \_\_\_\_\_

3) **Custom cell fixation requirements**

a. Cell fixant \_\_\_\_\_ (We typically use 10% buffered formalin.)

b. Temperature \_\_\_\_\_ (We typically use a controlled incubator (24°C +/-0.1°C; Select range from 24 - 60°C.)

c. Fixation time in hours \_\_\_\_\_ (We typically perform 24 hours incubation.)

**Turnaround: We usually deliver a FFPE Cell Block in 7 working days starting from fully grown and fixed cells.**

## Optional FFPE Services

1) **Need sections cut (i.e. scrolls) from the prepared block?** Yes or No

a. If yes, how many? \_\_\_\_\_

b. What section thickness \_\_\_\_  $\mu\text{m}$  (5 – 20  $\mu\text{m}$ , typically 20  $\mu\text{m}$ )

c. Please specify the position on the block to section \_\_\_\_\_ (e.g. top, middle or bottom of block)

Depending on thickness requested we can usually obtain 200 sections (5  $\mu\text{m}$  thick).

2) **Need slides prepared?** Yes or No If yes, how many? \_\_\_\_\_

3) **Need H&E staining of the slide(s)?** Yes or No If yes, how many? \_\_\_\_\_

4) **Need additional QC testing?** Yes or No If yes, please specify \_\_\_\_\_ (e.g. extract gDNA and quantification)

Call to inquire about **Custom Multiplex Cell Line Mixtures!**

[www.ONCOREF.com](http://www.ONCOREF.com)